

Super Pfx DNA Polymerase

Single Cell Whole Genome Amplification Kit (MDA)

Catalog Number:

s665716 (24 rxns)

s665716 (96 rxns)

Storage conditions: Please send the kit on dry ice and store all components in a refrigerator at -20°C for up to 6 months immediately after receipt. For longer term storage, please store below -70°C.

Products content

Component	24 rxns	96rxns
SC-DNA Polymerase	48 µl	192 µl
SC-Reaction Buffer	1 ml	4 x 1 ml
Buffer D	1 ml	1.5 ml
Buffer N	1 ml	1.5 ml
DTT, 1 M	1 ml	1 ml
PBS	1 ml	1.5 ml

Products Introduction

The Single Cell Whole Genome Amplification Kit is an isothermal amplification system based on MDA, which can be used as a template for whole genome amplification of single cells or micro samples. The size of single-cell whole genome amplification products ranges from 2-100 kb, which can be widely used in second-generation sequencing, large-band copy number variation analysis, microsatellite analysis, qPCR analysis, gene chip analysis and so on. The Phi29 DNA polymerase used in this kit is a DNA polymerase cloned from phage, which has strong strand displacement activity and strand affinity, and can achieve continuous polymerization and extension of up to 100 kb in a single polymerization reaction, and its amplification products are suitable for a variety of downstream applications, and the Phi29 DNA polymerase has strong 3'-5' exonuclease activity, and can be widely used in second-generation sequencing, copy number variation analysis of large fragments, microsatellite analysis, qPCR analysis and gene chip analysis. Phi29 DNA polymerase also has strong 3'-5' exonuclease activity, which ensures high fidelity of DNA synthesis. Under normal conditions, a single reaction can produce more than 20 µg of genomic DNA with high coverage.

Bring your own instruments and reagents

centrifuges

Water bath or PCR instrument

Reaction Tubes: Low adsorption PCR tubes are recommended.


Gun Heads: High quality filtered gun heads are recommended to prevent contamination.

deionized water

caveat

1. The sensitivity of this product is very high, the experimental operation should be completed in the positive pressure of the ultra-clean bench, the concentration of the amplification reaction product is high, should be isolated to avoid aerosol contamination caused by the amplification product.
2. The use of low-quality samples as templates can affect the quality of the final amplification product, and the use of heavily degraded and fragmented DNA as starting samples should be avoided.

Operation flow diagram



procedure

Cells as templates for amplification

This protocol is suitable for whole genome non-discriminatory amplification using 1-1000 cells as a template. Freshly prepared cell samples should be used to ensure the integrity of the starting genome, and apoptotic cells should not be used.

1. This protocol is suitable for whole genome non-discriminatory amplification using 1-1000 cells as a template. Freshly prepared cell samples should be used to ensure the integrity of the starting genome, and apoptotic cells should not be used.

individual parts making up a compound	volumetric
Buffer D	33 μl
DTT, 1 M	3 μl
total volume	36 μl

2. Add 4 μl of cell sample (resuspended in PBS) to the PCR tube. If the sample volume is less than 4 μl, make up to 4 μl using PBS.
3. Add 3 μl of Buffer D2, flick the wall of the tube to mix and centrifuge briefly to collect. Make sure that the cells are not adhering to the wall of the tube, and do not blow with the pipette to avoid the cell sample adhering to the pipette tip.
4. The samples were incubated at 65°C for 10 min.
5. Add 3 μl of Buffer N, flick the walls of the tube to mix and centrifuge briefly. Keep the sample on ice until ready for the next step.

6. Prepare the reaction mixture according to the table below, mix and centrifuge briefly.

individual parts making up a compound	volumetric
SC-Reaction Buffer	38 μ l
SC-DNA Polymerase	2 μ l
total volume	40 μ l

7. Immediately add 40 μ l of the reaction mixture to the prepared 10 μ l DNA sample (step 5), flick the walls of the tube to mix and collect by brief centrifugation.
8. Incubate at 30°C for 2 h. Incubation can be extended to increase yield if needed.
9. Incubate at 65°C for 5 min to inactivate SC-DNA Polymerase. Note: The amplified product is a high concentration of genomic DNA, please use water or TE to dilute to a suitable concentration for downstream experiments. The amplified product can be widely used in whole genome and exon sequencing, qPCR analysis, gene chip analysis, etc.

Genome as template amplification

This protocol is suitable for non-discriminatory genome-wide amplification of greater than 1 ng of purified genomic DNA as a template, but less starting DNA can be used if the genome is sufficiently complete and pure.

1. Prepare Buffer D1 and N1 (the volumes given in the table below are sufficient for 12 reactions, and may be stored at -20°C if not fully used in one experiment, but should not be stored for more than 3 months).

individual parts making up a compound	Buffer D1	Buffer N1
Buffer D	7 μ l	-
Buffer N	-	9 μ l
(of clothes) classifier for	25 μ l	51 μ l
number of washes	32 μ	60 μ l
total volume		

2. Add 2.5 μ l of DNA sample to the PCR tube. If the sample volume is less than 2.5 μ l, use water or TE to make up to 2.5 μ l.
3. Add 2.5 μ l of Buffer D1, flick the tube wall to mix and centrifuge briefly.
4. Incubate at room temperature (15-25°C) for 3 min.
5. Add 5 μ l of Buffer N1, flick the walls of the tube to mix and centrifuge briefly. Keep the sample on ice until ready for the next step.
6. Prepare the reaction mixture according to the table below, mix and centrifuge briefly.

individual parts making up a compound	volumetric
SC-Reaction Buffer	38 μ l
SC-DNA Polymerase	2 μ l
total volume	40 μ l

7. Immediately add 40 μ l of the reaction mixture to the prepared 10 μ l DNA sample (step 5), flick the walls of the tube to mix and collect by brief centrifugation.
8. Incubate at 30°C for 2 h. Incubation can be extended to increase yield if needed.
9. SC-DNA Polymerase was inactivated by incubation at 65°C for 5 min.

Note: The amplified product is a highly concentrated genomic DNA, please use water or TE to dilute it to a suitable concentration for downstream experiments. The amplified product can be

widely used in whole genome and exon sequencing, qPCR analysis, gene chip analysis and so on.

Usage

The following are examples of conventional PCR reaction systems and conditions, which should be improved and optimized according to the template, primer structure and fragment size.

1. PCR Reaction System All operations should be carried out on ice, and the components should be mixed well after thawing and stored at -20°C after use.

reagents	50 µL reaction system	final concentration
2×Super Pfx Buffer	25 µL	1×
dNTP Mix, 10 mM each	1.5-2.5 µL	300-500 µM each
Forward Primer, 10 µM	2 µL	0.4 µM
Reverse Primer, 10 µM	2 µL	0.4 µM
TemplateDNA, moderate amount	appropriate amount	<500 ng/50 µL
Super Pfx DNA Polymerase	0.5-0.75 µL	1-1.5 U/50 µL
ddH ₂ O	up to 50 µL	

2. PCR reaction conditions

move	temp	timing	
premutability	98°C	30 s-3 min	
denaturation	98°C	10-30 s	25-35 cycles
annealing (metallurgy)	According to the primer T _m	15-30 s	25-35 cycles
reach	72°C	3-5 kb/min	25-35 cycles
ultimate extension	72°C	5 min	

take note of

1) Priority is given to three-step amplification; if the three-step method fails to amplify the target product or if the primer T_m value is greater than 68° C, try the two-step method.

2) Denaturation: pre-denaturation of simple templates 98° C, 30s-1min, for complex templates, the pre-denaturation time can be extended to 3min.

3) Annealing: In general, the annealing temperature is 3-5°C lower than the T_m value of the primers. If the desired amplification efficiency cannot be obtained, the annealing temperature should be changed in a gradient to optimize the results; if a non-specific reaction occurs, the annealing temperature should be increased appropriately.

4) Extension: The extension time should be set according to the length of the amplified fragments and the complexity of the template, the amplification efficiency of this product is 3-5 kb/min, for long fragments and templates with high complexity it is recommended that 2-4kb/min.

5) Cycling times: the number of cycles can be set according to the downstream application of the amplification product. If the number of cycles is too small, the amplification amount will be insufficient, and if the number of cycles is too large, the chance of

mismatch will be increased, so the number of cycles should be minimized under the premise of guaranteeing the yield of the product.