

NGS TP Index Kit for Illumina (Index Primers Set II) Transposase Method for Second Generation Sequencing Multi-sample Primer Kit II (Illumina)

Catalog number: N665978 (384 rxns)

Storage condition: -20°C storage, dry ice transportation.

Products content

Component	384 rxns
Index N502-N522 Primers for Illumina	16 x 24 µl
Index N701-N729 Primers for Illumina	24 x 16 µl

Products Introduction

This kit is a companion to the transposase-based rapid DNA library construction kit, designed for Illumina platform library construction. It contains 16 N5 primers and 24 N7 primers, which can be used to prepare 384 different bipartite Index libraries. All reagents provided in the kit have been subjected to stringent quality control and functional validation to maximize the stability and reproducibility of library construction. The libraries can be used for sequencing on Illumina platforms such as HiSeq X-10/4000/2500/2000 and MiSeq.

Provide your own instruments, reagents and consumables

1. Magnetic frame: DynaMag^{TM-2} is recommended.
 2. DNA purification and recovery kit: It is recommended to use Kangwei DNA purification and recovery kit by magnetic bead method.
 3. DNA building kit: It is recommended to use the Kangwei Century transposase method second-generation sequencing rapid DNA building kit.
 4. Anhydrous ethanol.
 5. Reaction tubes: It is recommended to use low adsorption PCR tubes with 1.5 ml centrifuge tubes;
- Tip: It is recommended to use a high quality filter tip to prevent contamination of kits and library samples.

Pre-experiment Preparation and Important Notes

Please centrifuge briefly before opening the cap so that the liquid collects at the bottom of the tube to avoid cross-contamination between different primers.

procedure

For the use of the CombiVision Second Generation Sequencing Multisample Primer Kit, please follow the CombiVision Second Generation Sequencing Rapid DNA Library Kit protocol.

Index N502-N522 Primers for Illumina

	Index Primers for Illumina	Index
N502	5'-AATGATACGGCGACCACCGAGATCTA CACCTCTCTCTATTTCGTTCGGCAGCGTC-3'	CTCTCTAT
N503	5'-AATGATACGGCGACCACCGAGATCTA CACTATCCTCTTCGTTCGGCAGCGTC-3'	TATCCTCT
N505	5'-AATGATACGGCGACCACCGAGATCTA CACGTAAGGAGTCGTTCGGCAGCGTC-3'	GTAAGGAG
N506	5'-AATGATACGGCGACCACCGAGATCTA CACACTGCATATCGTTCGGCAGCGTC-3'	ACTGCATA
N507	5'-AATGATACGGCGACCACCGAGATCTA CACAAGGAGTATCGTTCGGCAGCGTC-3'	AAGGAGTA
N508	5'-AATGATACGGCGACCACCGAGATCTA CACCTAAGCCTTCGTTCGGCAGCGTC-3'	CTAAGCCT
N510	5'-AATGATACGGCGACCACCGAGATCTA CACCGTCTAATTCGTTCGGCAGCGTC-4'	CGTCTAAT
N511	5'-AATGATACGGCGACCACCGAGATCTA cactctctctccgctcgtcgcgcagcgtc-5'	TCTCTCCG

	Index Primers for Illumina	Index
N513	5'-AATGATACGGCGACCACCGAGATCTA cactcgactagtcgctcggcagcgtc-6'	TCGACTAG
N515	5'-AATGATACGGCGACCACCGAGATCTA CACTTCTAGCTTCGTTCGGCAGCGTC-7'	TTCTAGCT
N516	5'-AATGATACGGCGACCACCGAGATCTA CACCTAGAGTTCGTTCGGCAGCGTC-8'	CCTAGAGT
N517	5'-AATGATACGGCGACCACCGAGATCTA cacgcgtaagatcgtcgcgcagcgtc-3'	GCGTAAGA
N518	5'-AATGATACGGCGACCACCGAGATCTA CACCTATTAAGTCGTTCGGCAGCGTC-4'	CTATTAAG
N520	5'-AATGATACGGCGACCACCGAGATCTA CACAAGGCTATTCGTTCGGCAGCGTC-5'	AAGGCTAT
N521	5'-AATGATACGGCGACCACCGAGATCTA cacgagccttatcgtcgtcggcgtc-6'	GAGCCTTA
N522	5'-AATGATACGGCGACCACCGAGATCTA CACTTATGCGATCGTTCGGCAGCGTC-7'	TTATGCGA

Index N701-N729 Primers for Illumina

	Index Primers for Illumina	Index
N701	5'-CAAGCAGAAGACGGCATAACGAGAT TCGCCTTAGTCTCGTGGGCTCGG-3'	TAAGGCGA
N702	5'-CAAGCAGAAGACGGCATAACGAGAT CTAGTACGGTCTCGTGGGCTCGG-3'	CGTACTAG
N703	5'-CAAGCAGAAGACGGCATAACGAGAT TTCTGCCTGTCTCGTGGGCTCGG-3'	AGGCAGAA
N704	5'-CAAGCAGAAGACGGCATAACGAGAT GCTCAGGAGTCTCGTGGGCTCGG-3'	TCCTGAGC
N705	5'-CAAGCAGAAGACGGCATAACGAGAT AGGAGTCCGTCTCGTGGGCTCGG-3'	GGACTCCT
N706	5'-CAAGCAGAAGACGGCATAACGAGAT CATGCCTAGTCTCGTGGGCTCGG-3'	TAGGCATG
N707	5'-CAAGCAGAAGACGGCATAACGAGAT GTAGAGAGAGGTCTCGTGGGCTCGG-3'	CTCTCTAC
N710	5'-CAAGCAGAAGACGGCATAACGAGAT CAGCCTCGGTCTCGTGGGCTCGG-3'	CGAGGCTG
	Index Primers for Illumina	Index
N711	5'-CAAGCAGAAGACGGCATAACGAGAT TGCCTCTTGTCTCGTGGGCTCGG-3'	AAGAGGCA
N712	5'-CAAGCAGAAGACGGCATAACGAGAT tcctctacgtctcgtgggctcgg-3'	GTAGAGGA
N714	5'-CAAGCAGAAGACGGCATAACGAGAT TCATGAGCGTCTCGTGGGCTCGG-3'	GCTCATGA
N715	5'-CAAGCAGAAGACGGCATAACGAGAT cctgagatgtctcgtgggctcgg-3'	ATCTCAGG
N716	5'-CAAGCAGAAGACGGCATAACGAGAT TAGCGAGTGTCTCGTGGGCTCGG-3'	ACTCGCTA
N718	5'-CAAGCAGAAGACGGCATAACGAGAT GTAGCTCCGTCTCGTGGGCTCGG-3'	GGAGCTAC
N719	5'-CAAGCAGAAGACGGCATAACGAGAT TACTACGCGTCTCGTGGGCTCGG-3'	GCGTAGTA
N720	5'-CAAGCAGAAGACGGCATAACGAGAT AGGCTCCGTCTCGTGGGCTCGG-3'	CGGAGCCT

N721	5'-CAAGCAGAAGACGGCATAACGAGAT	TACGCTGC
N722	5'-CAAGCAGAAGACGGCATAACGAGAT CTGCGCATGTCTCGTGGGCTCGG-3'	ATGCGCAG
N723	5'-CAAGCAGAAGACGGCATAACGAGAT GAGCGCTAGTCTCGTGGGCTCGG-3'	TAGCGCTC
N724	5'-CAAGCAGAAGACGGCATAACGAGAT CGCTCAGTGTCTCGTGGGCTCGG-3'	ACTGAGCG
N726	5'-CAAGCAGAAGACGGCATAACGAGAT GTCTTAGGGTCTCGTGGGCTCGG-3'	CCTAAGAC
N727	5'-CAAGCAGAAGACGGCATAACGAGAT ACTGATCGGTCTCGTGGGCTCGG-3'	CGATCAGT
N728	5'-CAAGCAGAAGACGGCATAACGAGAT TAGCTGCAGTCTCGTGGGCTCGG-3'	TGCAGCTA
N729	5'-CAAGCAGAAGACGGCATAACGAGAT GACGTCGAGTCTCGTGGGCTCGG-3'	TCGACGTC