

# **Annealing Buffer for RNA Oligos (10X)**

# A743818

**Storage:** Store at -20 °C for up to 1 year.

#### Introduction:

Annealing Buffer for RNA Oligos (10X) is used for annealing complementary RNA oligos to form double-stranded RNA to synthesize siRNA after cell transfection. This product can ensure proper annealing by effective prevention of RNA oligos from forming hairpin structures. This product is easy-to-use. Users just need to mix RNA oligos and the Annealing Buffer for RNA Oligos (10X) in an appropriate ratio and place the reaction in a thermocycler to complete the annealing in 60 minutes.

This product is free RNase to avoid RNA degradation.

### **Usage method:**

- Dilute the RNA oligos to be annealed to 200μM with DEPC-treated water. Thaw the Annealing Buffer for RNA Oligos (10X) and mix well.
- 2. Add the reagents in order as follows and mix well.

Component	Volume
DEPC-treated Water	40μΙ
Annealing Buffer for RNA Oligos (10X)	10µl
RNA oligo A (200µM)	25µl
RNA oligo B (200µM)	25µl
Total volume	100μΙ

3. Set up the PCR instrument for annealing reaction.

Step	Tempe <mark>ra</mark> ture	Time
Denaturation	90℃	1 minute
Annealing	-0.1℃/ <mark>5s</mark> , to 25℃	60 minutes
Storage	<b>4</b> °C	ω

Note 1: The gradient temperature decreasing rate can also be set at -1°C/60s.

Note 2: When a thermocycler is not available, the PCR tube can be placed in a water bath with appropriated amount of 90°C water that can achieve the temperature decreasing from 90°C to 25°C within 1-2 hours naturally. However, the annealing efficacy will be inferior to that obtained using a thermocycler.

4. The annealed products can be directly used for cell transfection, or stored at -20  $^{\circ}$ C for future use.



### **Precautions:**

- 1. Annealing Buffer for RNA Oligos (10X) is for RNA oligos only, not for DNA oligos.
- 2. As RNA is very susceptible to degradation, extreme care must be taken to prevent RNase contamination during the course of experimentation.
- 3. This product is for R&D only. Not for drug, household, or other uses.
- 4. For your safety and health, please wear a lab coat and disposable gloves during the operation.

