

NP-40 Lysis Buffer (Contains inhibitors)

N748603

Storage: NP-40 Lysis Buffer store at 2-8°C,
Phosphatase Inhibitor Cocktail and PMSF Store at -20°C

Introduction:

NP-40 Lysis Buffer is a relatively mild cell and tissue lysis solution. The protein samples obtained by lysis with it can be used for routine Western blotting, immunoprecipitation (IP), and co-immunoprecipitation (co-IP), etc. For the protein samples obtained by lysis with NP-40 Lysis Buffer, the protein concentration can be determined using a BCA Protein Assay Kit. Due to the presence of a relatively high concentration of detergent, the Bradford method cannot be used to determine the protein concentration of the samples lysed by this lysis buffer.

Usage method:

1. Reagent Preparation

Take an appropriate amount of lysis buffer. Within a few minutes before use, add the protease inhibitor cocktail and/or phosphatase inhibitor cocktail to the lysis buffer at a ratio of 1:50. Alternatively, use 100 mM PMSF and adjust its final concentration to 1 mM.

2. Cell/Tissue Lysis

For adherent cells

- Remove the culture medium and wash the cells once with PBS, normal saline, or serum - free culture medium.
- Add the lysis buffer at a ratio of 150 - 250 μ L per well of a 6 - well plate. Pipette the lysis buffer up and down several times to ensure full contact between the lysis buffer and the cells. Usually, the cells will be lysed within 1 - 2 seconds after contact with the lysis buffer.
- If it is for ChIP (Chromatin Immunoprecipitation), continue the lysis on ice for another 10 minutes after the initial lysis.

For suspension cells

- Collect the cells by centrifugation and flick the tube forcefully with your finger to disperse the cells.
- Add the lysis buffer at a ratio of 150 - 250 μ L per well of a 6 - well plate equivalent of cells. Then gently flick the tube to fully lyse the cells.
- If the cell quantity is large, aliquot the cells into tubes at a density of 500,000 - 1,000,000 cells per tube before lysis. There should be no obvious cell pellet after full lysis.
- If it is for ChIP, continue the lysis on ice for another 10 minutes after the initial lysis.

For tissue samples

- Cut the tissue into small pieces with tissue scissors.

- b) Add the lysis buffer at a ratio of 150 - 250 μ L per 20 mg of tissue and homogenize the tissue using a glass homogenizer until fully lysed.
- c) If it is for ChIP, continue the lysis on ice for another 10 minutes after the initial lysis.

3. Protein Sample Collection

After full lysis, centrifuge the lysate at 10,000 - 14,000 g for 3 - 5 minutes. Collect the supernatant for subsequent operations such as PAGE (Polyacrylamide Gel Electrophoresis), Western blotting, and ChIP.

Precautions:

- 1. For the best usage effect, it is advisable to aliquot the reagent and store it at -20°C. Avoid excessive repeated freeze - thaw cycles.
- 2. All steps of lysing the samples should be carried out on ice or at 2 - 8°C.
- 3. For your own safety, please take protective measures such as wearing a lab coat and gloves before using the reagent.