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


Universal Protein Extraction Sample Kit effectively preserves protein post-translational modifications. (A-B) HeLa cells were serum starved for 16 h and stimulated with insulin (150 nM) for 5 min before protein extraction using the Total Protein Extraction Sample Kit. Total lysates (50 ug) were immunoblotted with anti phospho-mTOR (Ser2448) (A) and anti-mTOR (B) antibodies, respectively. (C) HeLa cells were lysed using the RIPA buffer or the Total Protein Extraction Sample Kit. Cell lysates were immunoblotted with anti-O-GlcNAc (glycosylation) antibody.

The Total Protein Extraction Sample Kit is effective in extracting proteins from tissues and cells. (A) HeLa cells and mouse brains (in duplicate) were lysed using the Total Protein Extraction Sample Kit. Total lysates (50 ug) were immunoblotted with anti-mTOR antibody. (B) The same lysates in (A) were immunoblotted with anti-GAPDH antibody. Note that the Total Protein Extraction Sample Kit is suitable for extracting proteins with a broad molecular weight range from both cells and tissues.

The Total Protein Extraction Sample Kit demonstrates superior performance in extracting large proteins. (A) HeLa cells were lysed using the Total Protein Extraction Sample Kit or RIPA buffer. Total lysates (50 ug) were immunoblotted with anti-mTOR antibody. (B) HT1080 cells were lysed using the Kit, processed for the indicated period, and immunoblotted with anti-SMRT antibody. Note that sonication fragmented SMRT to smaller peptides.