

RTL Reverse Transcriptase (Glycerol-free)

rp216167

Description:

RTL reverse transcriptase 2.0 is an RNA template-dependent DNA polymerase that lacks the 3' → 5' exonuclease activity and has RNase H activity. This enzyme can use RNA as a template to synthesize a complementary strand of DNA, which can be applied to first-strand cDNA synthesis, especially for RT-LAMP (loop-mediated isothermal amplification). Compared with RTL reverse transcriptase 1.0, the sensitivity is significantly improved, the thermal stability is stronger, and the reaction at 65°C is more stable. RTL reverse transcriptase 2.0 (glycerol free) can be used to prepare freeze-dried preparations, freeze-dried RT-LAMP reagents, etc

Storage and transportation

Storage at -25~-15°C and transportation < 0°C.

Unit definition

One unit incorporates 1 nmol of dTTP into acid-precipitable material in 20 minutes at 50°C using poly(A)•oligo(dT)25 as template-primer.

Composition of product

RTL Reverse Transcriptase 2.0 (Glycerol-free) (15U/μL)	0.1mL
10 × HH RTL Buffer	1.5mL
MgSO4(100 mM)	1.5mL

Quality Control Assays

Residual Activity of Endonuclease: A 50 μL reaction containing 1 μg of λDNA and 15 units of RTL2.0 incubated for 16 hours at 37°C shows same pattern as negative control by gel electrophoresis.

Residual Activity of Exonuclease: A 50 μL reaction containing 1 μg of Hind III digested λ DNA and 15 units of RTL2.0 incubated for 4 hours at 37°C shows same pattern as negative control by gel electrophoresis.

Residual Activity of Nickase: A 50 μL reaction containing 1 μg of supercoiled pBR322 and 15 units of RTL2.0 incubated for 4 hours at 37°C shows same pattern as negative control by gel electrophoresis.

Residual Activity of RNase: A 10 μL reaction containing 0.48 μg of MS2 RNA and 15 units of RTL2.0 incubated for 4 hours at 37°C shows same pattern as negative control by gel electrophoresis.

E.coli gDNA: Measured with E.coli specific HCD detection kits , 15 units of RTL2.0 contains less than 1 E. coli genome.

Reaction setup

Component	Volume
Template RNA ^a	optional
Oligo(dT)18~25 (50 μ M) or Random Primer Mix (60 μ M)	2 μ L
dNTP Mix (10 mM each)	1 μ L
RNase Inhibitor (40 U/ μ L)	0.5 μ L
RTL Reverse Transcriptase 2.0 (15 U/ μ L)	0.5 μ L
10 x HH RTL Buffer	2 μ L
Nuclease-free Water	Up to 20 μ L

Note: The recommended dosage of Total RNA is 1ng~1 μ g

The recommended dosage of mRNA was
50ng~100ng

Thermocycling Conditions for a routine reaction:

Temperature ($^{\circ}$ C)	Time
25 $^{\circ}$ C ^①	5 min
55 $^{\circ}$ C	10 min ^②
80 $^{\circ}$ C	10 min

①If Random Primer Mix is used, an incubation step at 25 $^{\circ}$ C.

②If target primer mix is used, an incubation step at 55 $^{\circ}$ C for 10~30 min.

RT-LAMP Protocol

Component	Volume	Fin. Con
Template RNA	optional	≥ 10 copies
dNTP Mix (10 mM)	3.5 μ L	1.4 mM
FIP/BIP Primers (25 $\times$)	1 μ L	1.6 μ M
F3/B3 Primers (25 $\times$)	1 μ L	0.2 μ M
LoopF/LoopB Primers (25 $\times$)	1 μ L	0.4 μ M
RNase Inhibitor (40 U/ μ L)	0.5 μ L	20U/Reaction
RTL Reverse Transcriptase 2.0 (15 U/ μ L)	0.5 μ L	7.5 U/Reaction
Bst V2 DNA Polymerase(8 U/ μ L)	1 μ L	8 U/Reaction
MgSO ₄ (100 mM)	1.5 μ L	6 mM (Total 8 mM)
10 x HH RTL Buffer(or 10 x HH Bst V2 Buffer*)	2.5 μ L	1 $\times$ (2 mM Mg ²⁺)
Nuclease-free Water	Up to 25 μ L	—

Mix by vortexing and centrifuge briefly to collect. Constant temperature incubation at 65 $^{\circ}$ C for 1 hour.

Note: The two buffers are interoperable and have the same composition

Precaution

- 1) This product will form a white solid when stored at -20 °C. Take it out from -20°C and put it on ice for about 10 minutes. After melting, it can be used by shaking and mixing.
- 2) The cDNA product could be stored at -20°C or -80°C or used immediately for PCR reaction.
- 3) To prevent RNase contamination, please keep the experimental area clean, and wear clean gloves and masks during operation.

