

β-Sodium Glycerophosphate Solution (5%)

Cat. No.: S1516055 | Pack size: 50 mL, 100 mL | Storage: 2-8°C, protected from light

Overview

β-Glycerophosphate Sodium Solution (5%) is mainly composed of β-glycerophosphate sodium and deionized water. It is one of the components of alkaline phosphatase staining solution and can also be used as an auxiliary component of protease inhibitors.

Protocol (for reference only)

Perform according to specific experimental requirements.

Specifications

Attribute	Value
Synonyms	β-Glycerol phosphate disodium salt solution (5%) β-Glycerophosphate Solution
Specifications & Purity	BioReagent, 5% in deionized water
Stability And Storage	Store at 2-8°C long term (12 months). Store in the dark.
Storage Conditions	Store at 2-8°C, Protected from light
Shipped In	Wet ice. This product requires cold chain shipping. Ground and other economy services are not available.

Contact & Global Offices

Whether you have a technical question, need help with a quotation, or want to inquire about an order, our regional teams are ready to assist. Please contact the office for your region; for general inquiries, the North American office is the corporate primary.

NORTH AMERICAN SALES, SUPPORT & GENERAL INQUIRIES Aladdin Scientific Corporation 14078 Meridian Parkway, Riverside, CA 92518, USA	EU SALES, LOGISTICS & LOCALIZED SUPPORT Aladdin Biochem Deutschland GmbH Westring 2, 33142 Büren, Germany
--	--



Phone: 1-833-552-7181

Sales: sales@aladdinsci.com

Customer Service: custserv@aladdinsci.com

Phone: +02951 9383958

Support: support.eu@aladdinsci.com

Limitations & Disclaimer

- For Research Use Only (RUO). Not for use in human or animal diagnostics, therapeutics, or in vivo applications. Not for food, cosmetic, or household use.
- This product is not a CE-marked in vitro diagnostic device under IVDR (EU) 2017/746 and is not an FDA-cleared device under 21 CFR. Use is restricted to verified businesses, institutions, and qualified professionals for research and development purposes.
- Where any kit component is classified as hazardous under CLP (EC) 1272/2008 or OSHA HCS (29 CFR 1910.1200), the product Safety Data Sheet (SDS) takes precedence over this document for handling, storage, and disposal information.
- Performance depends on sample type, sample condition, handling, and operator technique. Users are responsible for validating the product for their specific application.
- Aladdin product labels, SDS, COA, and approved specifications take precedence over this document. If product formulation, label, SDS, storage conditions, pack size, or quality specifications change, this document would be reviewed and reissued.