

Bovine Bile Salts

Cat. No.: B1522344 | Pack size: 100g; 250g | Storage: Room temperature, protected from light

Product Information

Item	Details
Available Pack Size(s)	100g; 250g
Specifications & Purity	BioReagent, for microscopy
Stability and Storage	Store at room temperature long term (36 months). Store in the dark.
Storage Conditions	Protected from light, Room temperature
Shipped In	Normal

Product Description

Bovine Bile Salts are sodium bile salts refined from bovine bile. They act by reducing the surface tension at the interface between the culture medium and bacterial pellicle, disrupting cell membrane structure and accelerating bacterial autolysis. They are widely used as common selective inhibitors in microbial culture media to inhibit the growth of Gram-positive bacteria.

Performance Indicators

Test Item	Standard Requirement
Appearance	White powder with a pungent bitter taste.
Clarity	Aqueous solution is clear and free of foreign impurities.
pH	7.0-9.0
Cholic Acid Content	≥60.0%
Residue on Ignition	≤17.0%
Moisture	≤5.0%

Precautions

1. Take precautions against dust during weighing. Wear a mask to avoid respiratory discomfort.
2. Tighten the bottle cap immediately after use of the powder to prevent moisture absorption and caking. Store in a dark and dry place.
3. The shelf life of unopened products is three years. The shelf life of opened products varies depending on storage conditions.

Contact & Global Offices

Whether you have a technical question, need help with a quotation, or want to inquire about an order, our regional teams are ready to assist. Please contact the office for your region; for general inquiries, the North American office is the corporate primary.

NORTH AMERICAN SALES, SUPPORT & GENERAL INQUIRIES

Aladdin Scientific Corporation

14078 Meridian Parkway, Riverside, CA 92518,
USA

Phone: 1-833-552-7181

Sales: sales@aladdinsci.com

Customer Service: custserv@aladdinsci.com

EU SALES, LOGISTICS & LOCALIZED SUPPORT

Aladdin Biochem Deutschland GmbH

Westring 2, 33142 Büren, Germany

Phone: +02951 9383958

Support: support.eu@aladdin.com

Limitations & Disclaimer

- For Research Use Only (RUO). Not for use in human or animal diagnostics, therapeutics, or in vivo applications. Not for food, cosmetic, or household use.
- This product is not a CE-marked in vitro diagnostic device under IVDR (EU) 2017/746 and is not an FDA-cleared device under 21 CFR. Use is restricted to verified businesses, institutions, and qualified professionals for research and development purposes.
- Where any product component is classified as hazardous under CLP (EC) 1272/2008 or OSHA HCS (29 CFR 1910.1200), the product Safety Data Sheet (SDS) takes precedence over this document for handling, storage, and disposal information.
- Performance depends on sample type, sample condition, handling, and operator technique. Users are responsible for validating this product for their specific application.



ALADDIN SCIENTIFIC

Reagents • Kits • Standards

Product Instructions
For Research Use Only
www.aladdinsci.com

- Aladdin product labels, SDS, COA, and approved specifications take precedence over this document. If product formulation, label, SDS, storage conditions, pack size, or quality specifications change, this document would be reviewed and reissued.