

Protease Inhibitor Cocktail

P665818-1mL

Store at -20°C, avoid freeze/thaw cycle.

Contents:

Cat No.	Component	Size	Storage
P665818-1mL	Protease Inhibitor Cocktail (100×)	1mL	-20°C

Notes:

1. Wear proper personal protective equipment (lab coats and disposable gloves) when working with this protease inhibitor cocktail.
2. Protease inhibitors harm the mucous membranes of the respiratory tract, eyes, and skin; inhalation, ingestion, or dermal absorption poses a lethal risk.

Operation steps:

1. Equilibrate the bottle to room temperature.
2. Vortex the tube before use to ensure a homogeneous suspension. Puncture the foil seal with a pipette tip.
3. Immediately before use, add 10µL of Protease Inhibitor Cocktail per milliliter of sample directly to the cell lysis buffer or extract to produce a 1X final concentration.