
Protease Inhibitor Blend

Project number: P665818

Storage conditions: -20° C.

Product Content:

Component	P665818
	1 ml
Protease Cocktail (100×)	Inhibitor 1 ml

Caveats:

1. Protease inhibitors cause serious damage to respiratory mucous membranes, eyes, and skin, and can be fatal if inhaled, swallowed, or absorbed through the skin.
2. Wear a lab coat and disposable gloves when handling.

Operational Steps:

Take out the Protein Extraction Reagent for pre-cooling, and add Protease Inhibitor Cocktail at a ratio of 1:99 (e.g., 10 μ l Protease Inhibitor Cocktail to 990 μ l of Extraction Reagent) to make a 1× working solution in the Extraction Reagent.