

Universal Virus Sample DNA/RNA Preservation Solution

Cat. No.: V751611 | Pack size: 100 mL; 500 mL | Storage: Room temperature

Product Information

Item	Description
Specifications & Purity	BioReagent, ready-to-use, sterile, RNase free
Stability And Storage	Store at Room temperature for 24 months. (under dry conditions)
Storage Conditions	Room temperature
Shipped In	Normal

Product Description

The preservation solution adopts an efficient cracking formula, which can quickly release and inactivate virus particles, make the sample lose infectivity and prevent the degradation of virus nucleic acid. The virus DNA/RNA was stored at room temperature for more than one month, which solved the problem of mold pollution after sampling.

Sample saving step

1. Oral swab and nasopharyngeal swab: 3ml preservation solution is recommended for a sampling swab;
2. Tissue homogenate: it is suggested to add preservation solution according to 1:5 (that is, 1ml tissue homogenate is added with 5ml preservation solution);
3. Saliva, whole blood, serum and plasma: it is recommended to add preservation solution in a ratio of 1:1 (that is, 1ml of saliva, whole blood, serum and plasma is added with 1ml of preservation solution);
4. Feces: It is recommended to add preservation solution according to 1:7 (that is, 1g of feces is added with 7ml of preservation solution).

Contact & Global Offices

Whether you have a technical question, need help with a quotation, or want to inquire about an order, our regional teams are ready to assist. Please contact the office for your region; for general inquiries, the North American office is the corporate primary.

<p>NORTH AMERICAN SALES, SUPPORT & GENERAL INQUIRIES</p> <p>Aladdin Scientific Corporation</p> <p>14078 Meridian Parkway, Riverside, CA 92518, USA</p> <p>Phone: 1-833-552-7181</p> <p>Sales: sales@aladdinsci.com</p> <p>Customer Service: custserv@aladdinsci.com</p>	<p>EU SALES, LOGISTICS & LOCALIZED SUPPORT</p> <p>Aladdin Biochem Deutschland GmbH</p> <p>Westring 2, 33142 Büren, Germany</p> <p>Phone: +02951 9383958</p> <p>Support: support.eu@aladdinsci.com</p>
---	--

Limitations & Disclaimer

- For Research Use Only (RUO). Not for use in human or animal diagnostics, therapeutics, or in vivo applications. Not for food, cosmetic, or household use.
- This product is not a CE-marked in vitro diagnostic device under IVDR (EU) 2017/746 and is not an FDA-cleared device under 21 CFR. Use is restricted to verified businesses, institutions, and qualified professionals for research and development purposes.
- Where any kit component is classified as hazardous under CLP (EC) 1272/2008 or OSHA HCS (29 CFR 1910.1200), the product Safety Data Sheet (SDS) takes precedence over this document for handling, storage, and disposal information.
- Performance depends on sample type, sample condition, handling, and operator technique. Users are responsible for validating the kit for their specific application.
- Aladdin product labels, SDS, COA, and approved specifications take precedence over this document. If product formulation, label, SDS, storage conditions, pack size, or quality specifications change, this document would be reviewed and reissued.