

Lactophenol Blue Solution

M767027

Store at room temperature in the dark.

Introduction:

Lactophenol blue solution can be used for the detection of fungal staining. It can stain fungi dark blue.

Usage method:

1. Drop a drop of lactophenol blue solution on a glass slide.
2. Use an inoculating needle to carefully spread the fungal culture medium evenly in the staining solution.
3. Slowly cover the cover slip from the edge of the liquid. Avoid the formation of air bubbles under the cover slip and let it stand for 5 minutes.
4. If necessary, seal the four sides of the cover slip with nail polish.
5. Observe under a microscope.

Staining result:

Fungi are stained dark blue.

Matters needing attention:

1. It is recommended to take 1-2 samples for a preliminary experiment when using this product for the first time.
2. For your own safety, please take precautions such as wearing a lab coat and gloves before using the reagent.