



ELISA Stop Solution (Ready-to-use)

E1372412

Storage room temperature

Introduction

ELISA Stop Solution is a proprietary solution used to terminate the peroxidase/TMB reaction for ELISA applications. The TMB substrate reacts with immobilized horseradish peroxidase (HRP) conjugated secondary antibodies to produce a blue solution. Color intensity is an indication of analyte level. After attaining the desired intensity, the reaction is terminated by addition of ELISA Stop Solution. Upon addition of ELISA Stop Solution on the color turns from blue to yellow.

Instructions for Use

Following TMB incubation, add an amount of ELISA Stop Solution equal to the amount of TMB substrate added to each well of the ELISA plate (typically 100 μ L per well). Use at full strength.