

---

## Chymotrypsin from porcine pancreas

### C1501322

---

#### Chymotrypsin

Source: Porcine pancreas

Chymotrypsin (chymotrypsin), endopeptidase. White or light yellow crystalline or amorphous powder, soluble in water, insoluble in organic solvents. When applied to protein, selective hydrolysis L-tyrosine and the carboxyl group of L-phenylalanine peptide bond formation.

#### Application

As pharmaceutical raw materials manufacturing injection. it always used to general trauma, wound healing that after surgery, anti-inflammatory, prevent local edema, hematocele, sprain hematoma, partial swelling that after breast surgery, tympanitis, rhinitis and soon. Use for the removal of cataract, Loose the lashes muscle ligament, reduce the cyst rupture and the hurt of retina. Besides, the combination use with the Trypsin and used as a broad-spectrum anti-inflammatory drugs.

As an aid in the medication of antibiotics and antibacterial drug, widely used in the field of animal medicine and animal health..

#### Storage Conditions

Store at 2-8°C long term (1 year). Store in the dark.