

Tannic Acid Aqueous Solution (5%)

T1516054

Storage: Room Temperature.

Shipping: Normal.

Introduction:

Tannic acid, also known as tannin, is widely used in tanning preparations. Its CAS number is 1401-55-4 and its molecular weight is 1701.20.

Tannic Acid Aqueous Solution (5%) consists of 5% tannic acid and deionized water. This reagent is for research use only. Not for clinical diagnosis or other purposes.

Protocol (for reference only):

Operate according to the specific experimental requirements.

Precautions:

1. All reagents mentioned above are irritating. Please handle with care.
2. For your safety and health, wear laboratory coat and disposable gloves during operation.