
DNA Phenol Reagent

Cat. No.: D1518314 | Pack size: 250 mL | Storage: Store at 2-8°C; protect from light

Overview

This reagent is a DNA extraction phenol reagent saturated with Tris-HCl (pH 8.0). It appears as a light yellow, clear liquid with a pungent odor. The upper layer consists of Tris-HCl buffer with pH > 7.8, while the lower layer is the phenolic phase, which contains 8-hydroxyquinoline as an antioxidant. This product is suitable for DNA isolation.

This reagent is for research use only and is not intended for clinical diagnosis or other purposes.

Materials to Prepare

- Tissue Lysis Buffer.
- Tris-HCl (1 M, pH 8.0).
- Chloroform.
- Isoamyl alcohol.
- Isopropanol.
- 75% Ethanol.
- TE Buffer.

Procedure (For Reference Only)

Take an appropriate amount of Tissue Lysis Buffer and follow the steps below for DNA extraction:

1. Adjust pH: Ensure the pH of the lysis buffer is around 8.0 (add 0.1 volume of 1 mol/L Tris-HCl, pH 8.0 for adjustment if necessary).
2. Organic Extraction: Sequentially add the following reagents:
 - 0.5 volume of DNA Phenol Reagent (Note: Allow the reagent to stand before use; use the lower organic phase. The upper layer is the Tris-HCl solution, which prevents oxidation of the lower organic phase).
 - 0.5 volume of Chloroform:Isoamyl alcohol (24:1)

Vortex for 10 seconds to mix thoroughly, then incubate on ice for 15 minutes.

3. Centrifugation and Separation: After centrifugation, carefully transfer the upper aqueous phase to a new centrifuge tube.

4. DNA Precipitation: Add an equal volume of isopropanol, mix well, and precipitate the DNA.

5. Washing and Dissolution: Wash the precipitate once with 75% ethanol, air-dry briefly, and dissolve in an appropriate volume of TE buffer for subsequent use.

Precautions

1. DNA Phenol Reagent is highly corrosive. Avoid skin contact or inhalation. If the reagent turns red or brown, it indicates oxidation and should not be used.
2. For your safety and health, please wear a lab coat and disposable gloves during operation.
3. Use the reagent as soon as possible after opening to prevent affecting subsequent experimental results.

Specifications

Attribute	Value
Specifications & Purity	BioReagent, Suitable for molecular biology
Stability And Storage	Store at 2-8°C long term (12 months). Store in the dark.
Storage Conditions	Store at 2-8°C; protect from light
Shipped In	Wet ice. This product requires cold chain shipping. Ground and other economy services are not available.

Contact & Global Offices

Whether you have a technical question, need help with a quotation, or want to inquire about an order, our regional teams are ready to assist. Please contact the office for your region; for general inquiries, the North American office is the corporate primary.

NORTH AMERICAN SALES, SUPPORT & GENERAL INQUIRIES

Aladdin Scientific Corporation

14078 Meridian Parkway, Riverside, CA 92518, USA

Phone: 1-833-552-7181**Sales:** sales@aladdinsci.com**Customer Service:** custserv@aladdinsci.com

EU SALES, LOGISTICS & LOCALIZED SUPPORT

Aladdin Biochem Deutschland GmbH

Westring 2, 33142 Büren, Germany

Phone: +02951 9383958**Support:** support.eu@aladdinsci.com

Limitations & Disclaimer

- For Research Use Only (RUO). Not for use in human or animal diagnostics, therapeutics, or in vivo applications. Not for food, cosmetic, or household use.
- This product is not a CE-marked in vitro diagnostic device under IVDR (EU) 2017/746 and is not an FDA-cleared device under 21 CFR. Use is restricted to verified businesses, institutions, and qualified professionals for research and development purposes.
- Where any kit component is classified as hazardous under CLP (EC) 1272/2008 or OSHA HCS (29 CFR 1910.1200), the product Safety Data Sheet (SDS) takes precedence over this document for handling, storage, and disposal information.
- Performance depends on sample type, sample condition, handling, and operator technique. Users are responsible for validating the product for their specific application.
- Aladdin product labels, SDS, COA, and approved specifications take precedence over this document. If product formulation, label, SDS, storage condition, pack size, or quality specifications change, this document would be reviewed and reissued.