



15 mg/ml D-Luciferin, Free Acid Stock Solution

Instructions

1. For a 15 mg/ml solution, begin by adding 233.36 mg potassium bicarbonate or 195.82 mg sodium bicarbonate to 66.6 ml sterile water.
2. Add 1.0 grams of D-Luciferin (Aladdin Catalog # [D115508](#)) into solution.
3. Gently mix by inversion until fully dissolved.
4. Aliquot into tubes.
5. Freeze solution at -80°C (Good for up to 1 year).

