

Rat Serum

R501133

Storage: Store at -20°C and avoid repeated freezing and thawing.

Introduction:

Rat serum is processed using blood from healthy rats as raw material, prepared through sterile collection, separation, microporous filtration, and aliquoting.

Usage Instructions:

1. Use according to specific experimental requirements.

Precautions:

1. Avoid repeated freezing and thawing. Repeated freezing and thawing may cause precipitation, which can be removed by centrifugation at 3500 rpm for 10 minutes under sterile conditions.
2. Pay attention to sterile operation.