



Rheumatoid Factor IgM from Human Plasma

np001022

Store at 2-8°C long term (>1 year). Upon receipt, it is recommended to aliquot.

Product Datasheet

Form: Liquid.

Buffer: Liquid in 20 mM sodium phosphate, pH 7.4, with 150 mM NaCl and 0.09% NaN₃.

Molecular Weight: 970,000 Da.

Protein Determination: Extinction Coefficient (E) = 1.18 (0.1% at 280 nm, 1 cm pathway).

Source: Human Plasma.

Testing: Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request.