


Gel Scan of Representative Lot

Lipoproteins, High Density, Human Plasma, Suitable for Cell Culture

No.np001038



SDS-PAGE:

12% BIO-RAD Tris HCl gel

1. LDL - 10 µg(Reduced/heat)
2. LMW Standard

Molecular Weight: 175,000-360,000

Considerations:

Storage Conditions: 2-8°C, away from light, Do Not Freeze

Stability: Open under a laminar flow hood; <= 7 days once opened

Buffer: 150 mM NaCl, pH 7.4, and 0.01% EDTA

Protein Determination:

Modified Lowry Based on BSA as standard

Density: 1.063-1.210 g/ml;

Composition: 55-45% lipid, 45-55% Protein

Sterility Testing:

Using Turbidity Growth/No Growth as the Indicator Thioglycollate Broth (THIO) for yeast and mold growth Tryptic Soy Broth (TSB) for bacterial growth

Endotoxin Level: <= 10 EU/mg by Pyrosate LAL clot test

Physical Specifications:

Liquid;

Purity: Single arc by IEP against antisera to whole human serum. Essentially free of other plasma lipoproteins as determined by electrophoresis using a SPIFE Vis Cholesterol gel kit for lipids and Coomassie Blue for proteins. Greater than 95% of total lipoprotein content by electrophoresis.

Product Datasheet

Product: Lipoproteins, High Density, Human Plasma, Suitable for Cell Culture

Number: np001038

Description: Liquid

Buffer: 150 mM NaCl, pH 7.4, and 0.01% EDTA

Purity: Single arc by IEP against antisera to whole human serum. Essentially free of other plasma lipoproteins as determined by electrophoresis using a SPIFE Vis Cholesterol gel kit for lipids and Coomassie Blue for proteins. Greater than 95% of total lipoprotein content by electrophoresis.

Sterility Testing: Using Turbidity Growth/No Growth as the Indicator Thioglycollate Broth (THIO) for yeast and mold growth Tryptic Soy Broth (TSB) for bacterial growth

Molecular Weight: MW 175,000-360,000

Density: 1.063-1.210 g/ml

Composition: 55-45% lipid, 45-55% Protein

Endotoxin Level: <= 10 EU/mg by Pyrosate LAL clot test

Protein Determination: Modified Lowry Based on BSA as standard

Storage: 2-8°C, away from light, Do Not Freeze

Stability: Open under a laminar flow hood; <= 7 days once opened

Available Packaging: 1mg per vial; Larger aliquots upon request

Source: Human Plasma

Testing: Shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.

Certificate Of Origin

Product Name: Lipoproteins, High Density, Human Plasma, Suitable for Cell Culture

SaturatedCatalog Number: np001038

Lot Number: Information applies to all lots produced.

Manufacturing Information: The methods and controls used in this facility for manufacturing, processing, packaging and holding of our products are in accordance with ISO9001:2015 Quality Management System.

Raw Material Information: The product contains no materials directly of animal origin and has not been in contact with products of animal origin in the course of its purification.

Country of Origin (manufacturing): United States

Biological Source Material (COO): United States

Donor Consent: The raw material used in the purification of this product

has been obtained from an FDA licensed collection center. Further the donor or the donor's next of kin has given by signature full informed consent for the material to be used in commercial scientific research, including the sharing, transferring, or selling of samples to third parties involved in this research. The raw material was obtained in an ethical manner in full compliance with national legislation, and the donor identity is held as confidential information by the licensed collection center.