

Gel Scan of Representative Lot
Lipoprotein a, [Lp(a)], Human Plasma
No.np001040

Product Datasheet

Product: Lipoprotein a, [Lp(a)], Human Plasma

Number: np001040

Description: Liquid

Buffer: 1X PBS with 0.05% NaN₃

Purity: >= 95% of total lipoprotein content by electrophoresis. SDS-PAGE shows minor ApoA1 contaminant.

Molecular Weight: 280,000-700,000Da

Density: 1.063-1.10 mg/ml

Composition: 60-65% lipid, 27% protein, 8% carbohydrate

Storage: 2-8°C away from light. Do Not Freeze

Stability: >= 6 months remaining

Source: Human Plasma

Testing: Source material shown to be non-reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA-required tests.

Certificate Of Origin

Product Name: Lipoprotein a, [Lp(a)], Human Plasma

Number: np001040

Raw Material Information: The product contains no materials directly of animal origin and has not been in contact with products of animal origin in the course of its purification.

Country of Origin (manufacturing): United States

Biological Source Material (COO): United States