


Hemopexin, Rat Serum

np001085



SDS-PAGE:

4-12% Bis-Tris NuPAGE gel

1. Hemopexin Rat - 5 µg (Reduced/heat)
2. Hemopexin Rat - 10 µg (Reduced/heat)
3. Hemopexin Rat - 20 µg (Reduced/heat)
4. LMW Standard
5. Hemopexin Rat - 5 µg (Non-reduced/no heat)
6. Hemopexin Rat - 10 µg (Non-reduced/no heat)
7. Hemopexin Rat - 20 µg (Non-reduced/no heat)

Molecular Weight: 57,000

Protein Determination:

Extinction Coefficient (E) 0.1% at 280nm =1.97
1cm pathway

Physical Specifications:

Lyophilized

Purity: ≥95% by SDS-PAGE

Considerations:

Storage Conditions: ≤-20 °C

Stability > 1 year

Buffer: 11.9 mM phosphate, 137 mM NaCl, 2.7 mM KCl, pH 7.4

Product Datasheet

Product:Hemopexin, Rat Serum

Catalog Number:np001085

Description:Lyophilized

Buffer:11.9 mM phosphate, 137 mM NaCl, 2.7 mM KCl, pH 7.4

Purity: $\geq 95\%$ by SDS-PAGE.

Molecular Weight:57,000

Protein Determination:Extinction Coefficient (E) 0.1% at 280nm =1.97
1cm pathway

Storage: $\leq -20^{\circ}\text{C}$

Stability:> 1 year

Available Packaging:500 μg per vial - 1 mg per vial

Source:Prepared from Rat Serum