


Haptoglobin, Cynomolgus Monkey Plasma

Project number: np001094



SDS-PAGE:

4-12% Bis-Tris-HCl NuPAGE gel

1. Haptoglobin - 5 µg (non-reduced / no heat)
2. Haptoglobin - 10 µg (non-reduced / no heat)
3. Haptoglobin - 20 µg (non-reduced / no heat)
4. BMW Standard
5. Haptoglobin - 5 µg (reduced / heated)
6. Haptoglobin - 10 µg (reduced / heated)
7. Haptoglobin - 20 µg (reduced / heated)

Molecular Weight:

86,000

Protein Determination:

Extinction Coefficient (E) $\frac{0.1\% \text{ at } 280\text{nm}}{1\text{cm pathway}} = 1.20$

Physical Specifications:

Lyophilized

Purity: $\geq 95\%$ by SDS-PAGE

Considerations:

Storage Conditions: $\leq -20^\circ\text{C}$

Stability > 1 year

Buffer: Salt-free lyophilized solid

Product Datasheet

Product: Haptoglobin, Cynomolgus Monkey Plasma

Catalog Number: np001094

Description: Lyophilized

Buffer: Salt-free lyophilized solid

Purity: $\geq 95\%$ by SDS-PAGE

Molecular Weight: 86,000

Protein Determination: Extinction Coefficient (E) $0.1\% \text{ at } 280\text{nm} = 1.20$
1cm pathway

Storage: $\leq -20^{\circ}\text{C}$

Available Packaging: 1 mg per vial
Larger aliquots upon request

Source: Monkey Plasma