



---

## Chymotrypsin from Human Pancreas

**np001104**

---

**Storage:** Store at -20°C. Avoid freeze/thaw cycle.

### Product Description

Chymotrypsin, Human Pancreas is a pancreas secreted serine protease that hydrolyses proteins into peptides and amino acids depending on its active group Ser-195.

Increased values of this enzyme and/or its zymogen have been found in patients with cystic fibrosis.

Prepared from tissue shown to be non reactive for HbsAG, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA-required tests.

