

Anti-Elastase Polyclonal Ab produced in rabbit

Ab001146

Product Datasheet

Product: Anti-Elastase Polyclonal Ab produced in rabbit

Catalog No: Ab001146

Description: Liquid

Host: Rabbit

Form: Ig (Standard Antisera)

Buier: phosphate buffered saline with 0.05% NaN₃ as a preservative

Storage: 2-8°C

Protein Concentration: 5.1mg/mL

Working Dilution: > 1:16000 by indirect ELISA

Western Blot: At 1:1000 dilution reacts to 0.3 ug purified antigen when run under non-reducing conditions. Does not show a reaction when antigen is run under reducing conditions.

IEP: Single arc against human neutrophil granule extract and human neutrophil elastase.