

Anti-Cathepsin H Polyclonal Ab produced in rabbit

Ab001148

Product Datasheet

Product: Anti-Cathepsin H Polyclonal Ab produced in rabbit

Catalog No: Ab001148

Description: Liquid

Host: Rabbit

Form: Ig (Standard Antisera)

Buier: phosphate buffered saline with 0.05% NaN₃ as a preservative

Storage: 2-8°C

Protein Concentration: 21.43mg/mL

Working Dilution: > 1:512 by indirect ELISA

Western Blot: At 1:8000 dilution reacts with 0.06 ug and .15 ug purified human liver cathepsin H. Does not react with cathepsins B, D or L.