

Anti-Catalase Polyclonal Ab produced in rabbit

Ab001151

Product Datasheet

Product: Anti-Catalase Polyclonal Ab produced in rabbit

Catalog No: Ab001151

Description: Liquid

Host: Rabbit

Form: IgG

Buer: phosphate buered saline with 0.05% NaN₃ as a preservative

Storage: 2-8°C

Protein Concentration: 1.96mg/mL

Working Dilution: >1:2000 by indirect ELISA using ALPH-conjugated secondary antibody

Western Blot: At 1:16,000 dilution, antisera reacts to 0.75, 0.15, and 0.3 ug catalase and to 1 uL of red blood cell lysate using ALPH conjugated secondary antibody.

IEP: Single arc upon immunoelectrophoresis of red blood cell lysate.