

Alpha 1 Acid Glycoprotein, Cynomolgus Monkey Plasma

np001163

Molecular Weight:~55,000 Da by SDS-PAGE

Protein Determination:

Extinction Coefficient (E) $\frac{0.1\% \text{ at } 280\text{nm}}{1\text{cm pathway}} = 0.89$

Physical Specifications:Lyophilized

Purity: >= 95% by SDS-PAGE

Considerations:

Storage Conditions: <= -20°C or lower

Stability > 1 year

Buer: deionized water

Product Datasheet

Product:Alpha 1 Acid Glycoprotein, Cynomolgus Monkey Plasma

Catalog Number:np001163

Description:Lyophilized

Buer:deionized water

Purity:>=95% by SDS-PAGE

Molecular Weight:~55,000 Da by SDS-PAGE

Protein Determination:

Extinction Coefficient (E) $\frac{0.1\% \text{ at } 280\text{nm}}{1\text{cm pathway}} = 0.89$

Storage: <=-20°C

Stability:> 1 Year

Available Packaging:100µg per vial

Larger aliquots upon request

Source:Cynomolgus Monkey Plasma