

## dNTP Mix (10 mM each)

Item No. D665873 (5 mL)

Storage condition: -30 ~ -15°C

### Product content

individual parts making up a compound	D665873 5ml
dNTP Mix (10 mM each)	5 x 1 ml

### Product Introduction

This product is an equimolar mixture of dATP, dGTP, dCTP and dTTP (10mM each), which can be used with all heat-resistant polymerases. dNTP purity is over 98%.

### Usage

dNTP (10 mM each) can be used directly or diluted to the appropriate concentration with sterile ultrapure water adjusted to neutral (pH 7.5). It can be used for PCR reaction, Real-time PCR, DNA sequencing, molecular labeling, cDNA synthesis and so on.