



Universal Antibody Dilution Buffer

W666033

Storage: Store at 2-8°C for 12 months; Store at -20°C for 24 months.

Introduction:

Universal Antibody Dilution Buffer is suitable for antibody dilution in Western blot (WB), immunofluorescence (IF), and immunohistochemistry (IHC).

Antibodies diluted with this buffer can be stored at 4°C for at least one month and reused multiple times for WB, IF, and IHC. The diluted antibodies can be used as long as they remain clear. Minor precipitation can be removed by centrifuging. Discard the diluted antibody if significant turbidity or microbial growth occurs.

This product contains optimized amounts of BSA, detergents, and antibody stabilizers to minimize nonspecific binding and enhance the stability of diluted antibodies.

Precautions:

1. To ensure antibodies can be reused multiple times, it is recommended to store the diluted antibodies at 4°C to slow down their degradation.
2. This product is for research use only, not for drug, clinical diagnosis, food, or other uses.
3. For your safety and health, please wear a lab coat and disposable gloves during the operation.