

Primary Antibody Dilution Solution

A666079

Store at 2-8°C

Introduction:

Primary Antibody Dilution Buffer can be used for dilution and preparation of primary antibodies. Its main ingredients including goat serum, detergent and primary antibody stabilization reagent, have been optimized to reduce non-specific bindings of primary antibody and stabilize primary antibodies after dilution. The diluted primary antibodies in this product can be stored at 2-8°C for six months at least, and can be reused multiple times.

Matters needing attention:

1. To recycle antibodies for later use, we recommend incubating membranes in primary antibody solutions at 4°C if possible to maintain the activity of primary antibodies.
2. For your safety and health, please wear a lab coat and disposable gloves during the operation.
3. This product is for R&D only. Not for drug, household, or other uses.