

NTP Mix (25mM each, Tris Buffered)

N748679

Store at -20°C

Introduction:

This NTP Mix (25mM each) is a mixed aqueous solution of ATP, CTP, GTP and UTP at 25mM each.

This product is free of DNase, RNase, phosphatase and protease, and can be used in different molecular biology and biochemistry applications, such as in vitro transcription with SP6, T3 and T7 RNA polymerases, RNA synthesis, and RNA amplification.

This product is prepared in 1mM Tris at pH 7.3-7.5.

Matters needing attention:

1. Can be thawed at room temperature. After thawing, it is advisable to store it in a cooler or in an ice bath. After use, it should be promptly stored at -20°C.
2. This product is intended for scientific research by professionals only and is not intended for clinical diagnosis or treatment. It must not be used for food or drugs and must not be stored in residential settings.
3. For your safety and health, please wear laboratory attire and disposable gloves when handling.