

## ddGTP, 10mM

**G759818**

---

Store at -20°C

### Introduction:

This product is a dideoxynucleotide, which acts as a chain-terminating agent for DNA polymerase. It is the primary substrate for nucleic acid mass spectrometry kits and can also be used in Sanger sequencing.

### Matters needing attention:

1. Can be thawed at room temperature. After thawing, it is advisable to store it in a cooler or in an ice bath. After use, it should be promptly stored at -20°C.
2. During the experimental process, please use DNase-free and RNase-free tips and centrifuge tubes.
3. This product is intended for scientific research by professionals only and is not intended for clinical diagnosis or treatment. It must not be used for food or drugs and must not be stored in residential settings.
4. For your safety and health, please wear laboratory attire and disposable gloves when handling.