

## **UTP (100mM, Tris Buffered)**

**U759858**

---

store at -20°C

### **Introduction:**

This UTP solution is supplied at a concentration of 100mM. It is free of DNase, RNase, phosphatase, and protease.

This product can be used in different molecular biology and biochemistry applications, such as in vitro transcription, RNA synthesis, and RNA amplification, etc.

This product is prepared in 1mM Tris and titrated to pH 7.3-7.5 with high-purity NaOH.

### **Matters needing attention:**

1. This product can be thawed at room temperature. Keep it on ice for use and store at -20°C immediately after use.
2. This product is for R&D only. Not for drug, household, or other uses.
3. For your safety and health, please wear a lab coat and disposable gloves during the operation.