



Catalog # B282503 Bis(triphenylphosphino)(2-methylphenyl)chloronickel(II), 99%




Technical Notes:

1. This complex is used as an air-stable nickel precatalyst for the amination of aryl chlorides, sulfamates, mesylates, and triflates.
2. This nickel precatalyst is used for the polymerization of halothiophenes.



Tech. Note (1)
Ref. (1)



Tech. Note (2)
Ref. (2)

References:

1. *Org. Lett.*, **2014**, 16, 220.
2. *Organometallics*, **2012**, 31, 2263.