




Catalog # D134292 Dichlorobis(triphenylphosphine)cobalt(II), 98%

Technical Notes:

1. Catalyst used for hydrostannations.
2. Catalyst used for alkyne-dihalomethaneamine couplings. – an efficient route for propargylamines.
3. Cobalt-catalyzed transformation of alkynyl C-H bond: aldehyde-alkyne-amine coupling.



References:

1. *Tetrahedron*, **2013**, 69, 4000.
2. *Tetrahedron Lett.*, **2012**, 53, 6199.
3. *SynLett.*, **2010**, 3, 475.