



Catalog #D196276 trans-Dichlorobis(triphenylphosphine)palladium(II)/potassium phosphate admixture
[Catkit single-use vials - 6.32 wt% Pd complex]

Note:

Weight-percent of components:

6.32 wt% palladium complex

93.68 wt% potassium phosphate

Technical Note:

1. Convenient, pre-weighed vial of palladium catalyst/base admixture useful for screening reactions. The vial contains 453mg of admixture, which will deliver 4 mole% of palladium catalyst and 2 equivalents of base, to a reaction using 1 mmole of substrate.