



Technical Notes:


1. Catalyst used for the addition of water to alkynes.
2. Highly efficient catalyst for the intermolecular hydroamination of alkynes.
3. Relay catalysis using a gold(I) complex/ Bronsted acid binary system for the synthesis of benzoxasiloles.
4. Gold(I) catalyzed intermolecular hydroarylation of alkenes with indoles.




Tech Note (1)
Ref. (1)




Tech Note (1)
Ref. (2)



Tech Note (2)
Ref. (3)



Tech Note (3)
Ref. (4)



Tech Note (4)
Ref. (5)

References:

1. *J. Am. Chem. Soc.*, **2003**, 125, 11925.
2. *Angew. Chem. Int. Ed.*, **2002**, 41, 4563.
3. *Org. Lett.*, **2003**, 5, 3349.
4. *RCS Adv.*, **2014**, 4, 6215.
5. *Chem. Eur. J.*, **2008**, 14, 8353.